



**THUASNE®**

---

THERAFIRM®

# **COMPRESSION SOLUTIONS**

*Product Catalog*




## WINGS FOR YOUR HEALTH

### Dear Customers,

From individuals seeking a little more comfort in their daily routines to clinicians and specialists dedicated to improving patient care, our mission is rooted in the belief that small things can create meaningful change. A simple pair of compression socks may seem ordinary, yet it can help someone move with greater ease, confidence, and freedom. Thuasne is committed to supporting people as they strive for better health and a better life. We work to empower wearers, providers, and caregivers alike. More than ever, we are honored to be your partner in delivering products that give patients wings for their health.

*Your Thuasne® North America Team*

A large industrial manufacturing facility, likely a textile mill, featuring rows of circular knitting machines. The machines are white and blue, with various spools of yarn and components visible. The floor is polished and reflective, and the ceiling has industrial lighting fixtures. The overall scene is bright and organized, showcasing a modern manufacturing environment.

Thuasne owns a manufacturing facility in North Carolina and is committed to producing its circular knit products in the USA that allow a closer direct link to retailers and the ability to shorten our customers' supply chain management.

**Our plant has obtained ISO 13485 certification.**



## **Innovative and Experienced**

With over 175 years of experience and the know-how of over six generations of entrepreneurs, Thuasne is firmly committed to empower each patient and caregiver in their efforts to restore mobility and improve quality of life through non-invasive care. Thuasne continues to focus on supporting patients in the fields of orthopedics, medical compression, and women's health.

## **Quality Compression**

Thuasne is known throughout the industry for quality medical products, which include Therafirm compression. Our compression socks and hosiery are monitored and control tested to meet quality standards and regulations. A dedicated in-house Quality Assurance team verifies that each compression product is made to meet these rigorous standards. In addition, Thuasne's compression manufacturing facility has been ISO 13485 certified since 2012. What this means is that the compression socks we provide are designed and produced for use in the prevention and treatment of diseases or other medical conditions. These medical devices help improve lives.

## **Our Brands**

**Ease**

**12**

**Core-Spun**

**26**

**Core-Sport**

**32**

**TheraSport**

**34**

**Therafirm**

**36**

**Compression  
Solutions**

**44**

# TABLE OF CONTENTS

- 6 What is Compression Therapy?**
- 6 Disclaimers Key**
- 7 Compression Levels Chart**
- 8 Product Selection Matrix**
- 10 True Gradient Compression from Thuasne®**
- 10 Care Instructions**
- 11 Measuring for Compression**
- 11 Donning Instructions**
- 12 Ease Gradient Compression Hosiery**
- 14 Size Charts**
- 15 Ease Opaque Knee Highs**
  - 15–20mmHg
  - 20–30mmHg
  - 30–40mmHg
- 15 Ease Opaque Knee Highs with Silicone Dot Band**
  - 20–30mmHg
  - 30–40mmHg
- 16 Ease Opaque Open-Toe Knee Highs**
  - 15–20mmHg
  - 20–30mmHg
  - 30–40mmHg
- 16 Ease Opaque Trouser Socks**
  - 15–20mmHg
  - 20–30mmHg
  - 30–40mmHg
- 17 Ease Opaque Thigh Highs**
  - 15–20mmHg
  - 20–30mmHg
  - 30–40mmHg
- 18 Ease Opaque Open-Toe Thigh Highs**
  - 15–20mmHg
  - 20–30mmHg
  - 30–40mmHg
- 18 Ease Opaque Pantyhose**
  - 15–20mmHg
  - 20–30mmHg
  - 30–40mmHg
- 19 Ease Opaque Open-Toe Waist Highs**
  - 15–20mmHg
  - 20–30mmHg
  - 30–40mmHg
- 20 Ease Sheer Knee Highs**
  - 15–20mmHg
  - 20–30mmHg
  - 30–40mmHg
- 20 Ease Sheer Open-Toe Knee Highs**
  - 15–20mmHg
  - 20–30mmHg
- 21 Ease Sheer Thigh Highs**
  - 15–20mmHg
  - 20–30mmHg
  - 30–40mmHg
- 22 Ease Sheer Open-Toe Thigh Highs**
  - 15–20mmHg
  - 20–30mmHg
- 22 Ease Sheer Pantyhose**
  - 15–20mmHg
  - 20–30mmHg
  - 30–40mmHg
- 23 Ease Sheer Open-Toe Pantyhose**
  - 15–20mmHg
  - 20–30mmHg
- 24 Ease Bold Socks**
  - 15–20mmHg
  - 20–30mmHg
- 24 Ease Lymphedema Arm Sleeves**
  - 20–30mmHg
  - 30–40mmHg
- 25 Ease Lymphedema Gauntlets**
  - 20–30mmHg
  - 30–40mmHg
- 25 Ease Lymphedema Gloves**
  - 20–30mmHg

**26 Core-Spun Gradient Compression Socks**

**28 Size Charts**

**29 Core-Spun Crew Socks**

10-15mmHg

**29 Core-Spun Cushioned Knee High Socks**

15-20mmHg

20-30mmHg

**30 Core-Spun Knee High Socks**

10-15mmHg

15-20mmHg

20-30mmHg

30-40mmHg

**31 Core-Spun Thigh High Socks**

15-20mmHg

20-30mmHg

30-40mmHg

**32 Core-Sport® Athletic Gradient Compression Socks**

**32 Size Chart**

**33 Core-Sport Athletic Socks**

15-20mmHg

20-30mmHg

**34 TheraSport® Athletic Gradient Compression Socks + Leg Sleeves**

**34 Size Chart**

**35 TheraSport Athletic Socks**

15-20mmHg

20-30mmHg

**35 TheraSport Athletic Leg Sleeves**

15-20mmHg

20-30mmHg

**36 Therafirm® Gradient Compression Hosiery**

**38 Size Charts**

**40 Therafirm Knee Highs**

10-15mmHg

15-20mmHg

20-30mmHg

30-40mmHg

**40 Therafirm Open-Toe Knee Highs**

20-30mmHg

30-40mmHg

**41 Therafirm Trouser Socks**

10-15mmHg

15-20mmHg

20-30mmHg

30-40mmHg

**41 Therafirm Thigh Highs with Lace-Top**

15-20mmHg

20-30mmHg

**42 Therafirm Thigh Highs**

15-20mmHg

20-30mmHg

30-40mmHg

**42 Therafirm Open-Toe Thigh Highs**

20-30mmHg

30-40mmHg

**43 Therafirm Microfiber Footless Tights**

10-15mmHg

**43 Therafirm Pantyhose**

15-20mmHg

20-30mmHg

**44 Compression Solutions**

**44 Therafirm Ease Adjust Inelastic Compression Below Knee Wrap**

**45 Therafirm Anti-Embolism Stockings**

**46 Knee Interfaces**

**46 Compressogrip® Sleeves**

**47 Compressogrip® Tubular Elastic Bandages**

**48 Online Ordering System**

**49 Terms + Conditions**

## WHAT IS COMPRESSION THERAPY?


Compression allows patients to live better. Whether they need to wear compression to prevent DVT in travel, prevent tired, achy legs while they're working a long shift, or want to enhance endurance during a run, compression is the solution they seek. Compression therapy uses socks or stockings to apply external pressure to the limbs, improving circulation and reducing or preventing swelling.

### How Compression Works

True gradient compression garments are designed to deliver a controlled amount of pressure which is greatest at the ankle of the garment and gradually decreases toward the top of the stocking to provide maximum vein support.

### How Compression is Measured

Compression is measured in millimeters of mercury per square inch (mmHg), or how much pressure is required to elevate mercury. Higher levels of compression offer greater therapeutic values for more serious pathologies. Use the chart at the right for suggested usages.



### Disclaimers Key

\*The mean compression of an average ankle size.

^ Contraindications: Compression products should not be worn and are contraindicated if you have any of the following conditions: arterial insufficiency or degeneration, acute DVT (deep vein thrombosis), untreated congestive heart failure, untreated cancer, untreated infection, absent or severely impaired sensation, allergy to materials, severe or moderate peripheral arterial disease, severe cognitive impairment.

\*\*Contraindications: Severe arterial insufficiency, cutaneous infections, acute dermatitis, wet dermatosis, uncontrolled congestive heart failure, skin irritations, allergies to dyes.

† Contraindications: Compression products should not be worn and is contraindicated if you have any of the following conditions: Severe arterial insufficiency, cutaneous infections, acute dermatitis, wet dermatosis, uncontrolled congestive heart failure, skin irritations, allergies to dyes.

## Compression Levels Chart

<b>10–15mmHg*</b> Light Compression	<b>15–20mmHg*</b> Mild Compression	<b>20–30mmHg*</b> Moderate Compression	<b>30–40mmHg*</b> Firm Compression
Light support gradient compression is ideal for those individuals required to sit or stand for long periods of time. Designed to help improve circulation and prevent swelling.	Great for daily wear, travel, and sports, compression stockings help improve circulation without being too tight on your legs. 15–20mmHg compression garments come in a variety of materials to fit any lifestyle.	Ideal for varicose veins, spider veins, swelling (edema) and after surgery. You don't have to have a condition to wear 20–30mmHg. People looking for a bit more support than the 15–20mmHg usually choose to wear 20–30mmHg.	This stronger level of compression is recommended for blood clots (DVT), lymphedema, severe swelling (edema), severe varicose veins and post-surgery.
Improve circulation	Relief of heavy, tired, achy legs	Relief of heavy, tired, achy legs	Relief of heavy, tired, achy legs
Help prevent swelling	Mild varicose veins	Moderate varicose veins	Advanced varicose veins
Energize legs and feet	Poor blood circulation	Venous insufficiency	Severe chronic venous insufficiency
<p><b>What does mmHg stand for?</b></p> <p>The medical abbreviation mmHg stands for “millimeters of mercury”, which is a unit of measurement for pressure. It's the same unit of measurement for measuring blood pressure. Compression stockings use a range for how much pressure it provides, because the compression changes depending on your size.</p> <p><i>This is not a recommendation for a specific product or brands. This is general information and not intended to diagnose or treat patients. It is the provider's responsibility to ascertain the appropriateness of individual applications of compression therapy.</i></p> <p><b>Patients should always consult with their provider prior to wearing compression products.</b></p>	Mild to moderate swelling	Moderate swelling	Moderate swelling
	Mild edema	Minimal edema or lymphedema	Moderate edema or lymphedema
	DVT for long distance travelers	Prevention or management of DVT	Prevention, management or treatment of DVT
	Varicosities and swelling during pregnancy	Varicosities and swelling during pregnancy	Post-Thrombotic Syndrome
	Thrombosis prophylaxis	Superficial thrombophlebitis	Post-Phlebitis Syndrome
	Individuals subjected to immobility	Arterial insufficiency	Post vein stripping or sclerotherapy
	Individuals on their feet for extended periods of time	Post-op sclerotherapy or vein ligation	After healing of moderate ulcerations
	Leg discomfort	Post venous ulceration	

# PRODUCT SELECTION MATRIX

	Classifications	Ease Opaque	Ease Sheer	Ease Bold	Ease Lymphedema
10-15 mmHg	<p><b>C0</b> No visible or palpable signs of venous disease</p> <p><b>C1</b> Telangiectasies or Reticular Veins</p>				
15-20 mmHg	<p><b>C1</b> Telangiectasies or Reticular Veins</p> <p><b>C2</b> Varicose Veins Reticular veins <math>\geq</math> 3mm</p>	<p>Trouser Socks</p> <p>Knee Highs</p> <p>Thigh Highs</p> <p>Pantyhose/Waist Highs</p>	<p>Knee Highs</p> <p>Thigh Highs</p> <p>Pantyhose</p>	<p>Trouser Socks</p>	
20-30 mmHg	<p><b>C2</b> Varicose Veins Reticular veins <math>\geq</math> 3mm</p> <p><b>C3</b> Edema</p> <p><b>C4</b> Pigmentation or Eczema</p> <p><b>C5</b> Healed Venous Ulcer</p>	<p>Trouser Socks</p> <p>Knee Highs</p> <p>Thigh Highs</p> <p>Pantyhose/Waist Highs</p>	<p>Knee Highs</p> <p>Thigh Highs</p> <p>Pantyhose</p>	<p>Trouser Socks</p>	<p>Arm Sleeves</p> <p>Gauntlets</p> <p>Gloves</p>
30-40 mmHg	<p><b>C5</b> Healed Venous Ulcer</p> <p><b>C6</b> Active Venous Ulcer</p>	<p>Trouser Socks</p> <p>Knee Highs</p> <p>Thigh Highs</p> <p>Pantyhose/Waist Highs</p>	<p>Knee Highs</p> <p>Thigh Highs</p> <p>Pantyhose</p>		<p>Arm Sleeves</p> <p>Gauntlets</p>

**Core-Spun**

**Core-Sport**

**TheraSport**

**Therafirm**

Crew Socks  


Knee High Socks  


Knee High Socks  


Cushioned Knee High Socks  


Thigh High Socks  


Knee High Socks  


Cushioned Knee High Socks  



Thigh High Socks  


Knee High Socks  


Thigh High Socks  


Blank

Athletic Socks  



Athletic Socks  



Blank

Blank

Athletic Socks  


Athletic Leg Sleeves  


Athletic Socks  


Athletic Leg Sleeves  


Blank

Trouser Socks  


Knee Highs  


Microfiber Footless Tights  


Pantyhose  


Trouser Socks  


Knee Highs  


Thigh Highs  



Pantyhose  



Trouser Socks  



Knee Highs  



Thigh Highs  


Pantyhose  


Trouser Socks  


Knee Highs  


Thigh Highs  


Pantyhose  


# THERAFIRM®

**TRUE GRADIENT  
COMPRESSION  
FROM THUASNE**

## Comfortable

Our garments are made of soft fibers and offer comfortable support encouraging patients to be more compliant with their compression therapy. Moisture wicking fibers keep patients dry by evaporating moisture away from the skin. Soft, comfortable knit-in bands keep garments in place without causing uncomfortable indentations.

## Affordable

Compression affordability can limit or prevent compliance with prescribed compression therapy. High quality, medically effective Therafirm compression products can be offered at a lower price than other brands. Often the affordable price will encourage wearers to obtain several pairs to allow a normal laundering interval, and to provide a “wardrobe” for various occasions.

EASE  
BY THERAFIRM

core-spun  
BY THERAFIRM

core sport  
BY THERAFIRM

THERA SPORT  
BY THERAFIRM

THERAFIRM®

## Care Instructions

In order to maintain compression efficacy and keep garments looking like new, we suggest the following care instructions.

▶ **Hand wash**

Wash garments by hand in lukewarm or cold water with a mild detergent.

▶ **Do not bleach**

Bleach will degrade the fibers in your compression products and cause them to break down faster.

▶ **Line or drip dry and do not iron**

Prolonged exposure to excessive heat will more quickly break down the fibers in your compression products.

▶ **Do not wring**

Wringing the product can distort and wrinkle the fibers and affect the efficacy of your compression.

## Measuring for Compression


The key to a proper fit is in the measurement. Follow the below instructions and tips to ensure the size you select will be the most accurate size for your patient.

### Measurement Instructions

1. Measure the ankle circumference at the narrowest part of the ankle, above the ankle bone.
2. Measure the calf circumference at the fullest part of the calf.
3. For Knee Highs, measure the length from the floor to just below the bend of the knee.
4. For Thigh Highs and Pantyhose, measure the thigh circumference at the widest part of the thigh.

### Measurement Tips

- ▶ The best time for patients to be measured is when they first get up in the morning. Measurements taken later in the day are not as accurate due to swelling that occurs throughout the day.
- ▶ Measure the patient while he or she is standing. If the fitting cannot take place before the patient has been on his/her feet for some time, the patient's legs and feet should be elevated for 15 minutes prior to taking the measurements.
- ▶ Take snug measurements against the skin.
- ▶ Avoid measuring over socks or slacks.
- ▶ Never allow folds in stockings.
- ▶ Never remove any wound dressings.



## Donning Instructions

- Step 1. Insert hand and grab stocking at top of heel pocket.
- Step 2. While still holding heel pocket, turn top of stocking down toward toe.
- Step 3. Open stocking and slide foot in until toe and heel are positioned in place.
- Step 4. **Knee Highs:** Grasp top of stocking and pull up over ankle and calf. Position top approximately two fingers below the bend of knee.

**Thigh Highs:** Grasp top of stocking and pull up over ankle and calf. Position top of stocking at mid-thigh for Therafirm and top of leg for Ease.

**Pantyhose:** With crotch in a snug, comfortable position, stretch panty area high with both hands. Adjust waist band to normal position.


### Notes

1. Ease by Therafirm and Core-Spun by Therafirm products are very stretchy so it is important the garment is placed correctly (see Step 4). In addition, to ensure a proper fit and optimal comfort, adjust the garment so the material is evenly distributed on the leg. Check that there are no creases or wrinkles.
2. To avoid snags in hosiery, keep nails trimmed and refrain from wearing.
3. Donning gloves are recommended if patient needs additional help donning.
4. Stockings should not be too tight over dressings. Socks must fit the shape of the foot. Socks should have even compression or elasticity so that the sock takes on the anatomical shape of the foot without high compression (binding points) that may create discomfort or circulation obstruction. Observe the bare foot prior to donning the sock. After donning the sock, the foot shape should be consistent with that of the bare foot. **If pain, discomfort, or areas of binding occur, the stockings should be removed. Always consult with your physician before wearing gradient compression garments.**

### How to Measure Video Here



### How to Don Video Here



# EASE

BY THERAFIRM®

## GRADIENT COMPRESSION HOSIERY

The Ease by Therafirm premium gradient compression line is engineered using ultra-stretchy yarns for compression garments that are easy to put on and more comfortable to wear. Super soft, breathable material offers superior comfort for all-day wear. Moisture wicking fibers provide a comfortable coolness.

### Compression Levels Available

**15–20mmHg\***  
Mild Support

Prevents mild swelling

**20–30mmHg\***  
Moderate Support

Prevents mild to moderate swelling

**30–40mmHg\***  
Firm Support

Prevents moderate swelling

Refer to page 6–7 for more information on compression levels.

15–20mmHg is not available in Ease Lymphedema products.



Smooth band on open-toe products for comfort and durability



Thigh highs feature comfortable silicone band



Soft knit-in waistband to prevent rolling and enhance comfort

Lab-tested true gradient compression promotes better blood flow, prevents mild to moderate swelling, and relieves tired, achy legs and feet

Ultra-stretchy yarns for garments that are easy to put on and more comfortable to wear




Arm sleeves include a knit pattern at the inside of the elbow for a more precise fit to help the sleeve to stay in place, along with a soft, non-restrictive wrist

Soft, comfortable band is non-restrictive while staying in place all day

Opaque knee highs are also available with silicone band



Lab-tested true gradient compression gradually decreases toward the wrist to promote blood flow



Super soft, breathable material offers superior all-day comfort

Moisture wicking fibers provide dry and comfortable coolness



Smooth toe seam

Reinforced toe and heel for durability

Available in petite, short, and long lengths



Seamless, compressive fingers in the glove

Ease

Core-Spun

Core-Sport


TheraSport

Therafirm

Compression Solutions


# Ease Size Charts

## Socks + Knee Highs




Size	Ankle	Calf	Length
Small	6.5"-8.5" (17-22cm)	11"-16.5" (28-42cm)	Petite 13" or less (33cm)
Medium	8"-10" (21-26cm)	12"-17.5" (31-45cm)	Short 13"-16" (33-41cm)
Large	9"-11.5" (23-29cm)	13"-19" (33-49cm)	Long Over 16" (41cm)
X-Large	11"-15" (28-38cm)	17"-23" (43-58cm)	

## Arm Sleeves




Size	Wrist	Elbow	Axilla	Length
Small	5.5"-6.5" (14-16.5cm)	7.5"-9" (19-23cm)	8.5"-11.5" (21.5-29cm)	Regular 16" or less (41cm)
Medium	6.25"-7.25" (16-18.5cm)	9"-10.5" (23-26.5cm)	11"-14" (28-35.5cm)	Long Over 16" (41cm)
Large	7"-8" (18-20cm)	10.5"-12.25" (26.5-31cm)	13.5"-16.5" (34-42cm)	
X-Large	7.75"-9" (19.5-23cm)	12.25"-13.75" (31-35cm)	16"-19" (40.5-48cm)	

## Thigh Highs




Size	Ankle	Calf	Thigh	Length
Small	6.5"-8.5" (17-22cm)	11"-16.5" (28-42cm)	15.5"-21.5" (40-54cm)	Short 28.5" (73cm) or less
Medium	8"-10" (21-26cm)	12"-17.5" (31-45cm)	17.5"-25.5" (45-65cm)	Long Over 28.5" (73cm)
Large	9"-11.5" (23-29cm)	13"-19" (33-49cm)	21.5"-29.5" (55-75cm)	
X-Large	11"-15" (28-38cm)	17"-23" (43-58cm)	25.5"-33.5" (65-85cm)	

## Gauntlets + Gloves



Size	Palm Girth	Wrist Girth
Small	6.75"-7.5" (17-19cm)	5.5"-6.5" (14-16.5cm)
Medium	7.5"-8.25" (19-21cm)	6.25"-7.25" (16-18.5cm)
Large	8.25"-9" (21-23cm)	7"-8" (18-20cm)
X-Large	9"-9.75" (23-25cm)	7.75"-9" (19.5-23cm)

## Pantyhose + Waist Highs



Size	Ankle	Calf	Thigh	Hip	Length
Small	6.5"-8.5" (17-22cm)	11"-16.5" (28-42cm)	15.5"-21.5" (40-54cm)	to 42" (107cm)	Short 28.5" (73cm) or less
Medium	8"-10" (21-26cm)	12"-17.5" (31-45cm)	17.5"-25.5" (45-65cm)	to 48" (122cm)	Long Over 28.5" (73cm)
Large	9"-11.5" (23-29cm)	13"-19" (33-49cm)	21.5"-29.5" (55-75cm)	to 56" (143cm)	
X-Large	11"-15" (28-38cm)	17"-23" (43-58cm)	25.5"-33.5" (65-85cm)	to 64" (163cm)	



## Ease Opaque Knee Highs

Opaque knee highs feature a soft, knit-in band for comfortable all-day wear. A reinforced toe and heel provide durability.



Black Sand Natural Navy

### 15–20mmHg\* Mild Support\*\*


	Black	Sand	Natural	Navy
<b>Petite</b>	Small Petite	52065	52060	
	Medium Petite	52075	52070	
	Large Petite	52085	52080	
	X-Large Petite	52095	52090	
<b>Short</b>	Small Short	52015	52010	52013 52017
	Medium Short	52025	52020	52023 52027
	Large Short	52035	52030	52033 52037
	X-Large Short	52045	52040	
<b>Long</b>	Small Long	53015	53010	53013 53017
	Medium Long	53025	53020	53023 53027
	Large Long	53035	53030	53033 53037
	X-Large Long	53045	53040	

### 20–30mmHg\* Moderate Support\*\*

	Black	Sand	Natural	Navy
<b>Petite</b>	Small Petite	50065	50060	
	Medium Petite	50075	50070	
	Large Petite	50085	50080	
	X-Large Petite	50095	50090	
<b>Short</b>	Small Short	50015	50010	50013 50017
	Medium Short	50025	50020	50023 50027
	Large Short	50035	50030	50033 50037
	X-Large Short	50045	50040	
<b>Long</b>	Small Long	54015	54010	54013 54017
	Medium Long	54025	54020	54023 54027
	Large Long	54035	54030	54033 54037
	X-Large Long	54045	54040	


### 30–40mmHg\* Firm Support\*\*

	Black	Sand
<b>Petite</b>	Small Petite	56065 56060
	Medium Petite	56075 56070
	Large Petite	56085 56080
	X-Large Petite	56095 56090
<b>Short</b>	Small Short	56015 56010
	Medium Short	56025 56020
	Large Short	56035 56030
	X-Large Short	56045 56040
<b>Long</b>	Small Long	59015 59010
	Medium Long	59025 59020
	Large Long	59035 59030
	X-Large Long	59045 59040



## Ease Opaque Knee Highs with Silicone Dot Band

Opaque knee highs deliver soft, reliable comfort and feature a silicone-dotted top band for improved grip. A reinforced heel and toe add durability.



Black Sand

### 20–30mmHg\* Moderate Support\*\*

	Black	Sand
<b>Short</b>	Small Short	57115 57110
	Medium Short	57125 57120
	Large Short	57135 57130
	X-Large Short	57145 57140
<b>Long</b>	Small Long	57215 57210
	Medium Long	57225 57220
	Large Long	57235 57230
	X-Large Long	57245 57240

### 30–40mmHg\* Firm Support\*\*

	Black	Sand
<b>Short</b>	Small Short	58115 58110
	Medium Short	58125 58120
	Large Short	58135 58130
	X-Large Short	58145 58140
<b>Long</b>	Small Long	58215 58210
	Medium Long	58225 58220
	Large Long	58235 58230
	X-Large Long	58245 58240



## Ease Opaque Open-Toe Knee Highs

Open-toe knee highs are crafted for everyday wear, featuring a smooth, knit-in band for a comfortable fit that stays put. The open-toe design makes them ideal for open-toe footwear.



Black



Sand

### 15–20mmHg\* Mild Support\*\*

	Black	Sand
● Short	Small Short	51015 51010
	Medium Short	51025 51020
	Large Short	51035 51030
	X-Large Short	51045 51040
● Long	Small Long	51715 51710
	Medium Long	51725 51720
	Large Long	51735 51730
	X-Large Long	51745 51740

### 20–30mmHg\* Moderate Support\*\*

	Black	Sand
● short	Small Short	50615 50610
	Medium Short	50625 50620
	Large Short	50635 50630
	X-Large Short	50645 50640
● Long	Small Long	54615 54610
	Medium Long	54625 54620
	Large Long	54635 54630
	X-Large Long	54645 54640

### 30–40mmHg\* Firm Support\*\*

	Black	Sand
● Short	Small Short	56615 56610
	Medium Short	56625 56620
	Large Short	56635 56630
	X-Large Short	56645 56640
● Long	Small Long	59615 59610
	Medium Long	59625 59620
	Large Long	59635 59630
	X-Large Long	59645 59640



## Ease Opaque Trouser Socks

Versatile trouser socks deliver comfort in neutral, everyday colors. They offer a traditional rib pattern for men and a classic chevron pattern for women. A reinforced toe and heel hold up to everyday wear and tear.



Women's  
Chevron  
Black



Women's  
Chevron  
Khaki



Men's  
Ribbed  
Black



Men's  
Ribbed  
Khaki



Men's  
Ribbed  
Navy

See size charts on page 14.

**15-20mmHg\* Mild Support\*\***

	Women's Chevron Black	Women's Chevron Khaki	Men's Ribbed Black	Men's Ribbed Khaki	Men's Ribbed Navy	
<b>Short</b>	Small Short	51165	51169	52815	52819	52817
	Medium Short	51175	51179	52825	52829	52827
	Large Short	51185	51189	52835	52839	52837
	X-Large Short	51195	51199	52845	52849	52847
<b>Long</b>	Small Long	51265	51269	53815	53819	53817
	Medium Long	51275	51279	53825	53829	53827
	Large Long	51285	51289	53835	53839	53837
	X-Large Long	51295	51299	53845	53849	53847

**20-30mmHg\* Moderate Support\*\***

	Women's Chevron Black	Women's Chevron Khaki	Men's Ribbed Black	Men's Ribbed Khaki	Men's Ribbed Navy	
<b>Short</b>	Small Short	52165	52169	51415	51419	51417
	Medium Short	52175	52179	51425	51429	51427
	Large Short	52185	52189	51435	51439	51437
	X-Large Short	52195	52199	51445	51449	51447
<b>Long</b>	Small Long	52265	52269	55415	55419	55417
	Medium Long	52275	52279	55425	55429	55427
	Large Long	52285	52289	55435	55439	55437
	X-Large Long	52295	52299	55445	55449	55447

**30-40mmHg\* Firm Support\*\***

	Men's Ribbed Black	Men's Ribbed Khaki	Men's Ribbed Navy	
<b>Short</b>	Small Short	57415	57419	57417
	Medium Short	57425	57429	57427
	Large Short	57435	57439	57437
	X-Large Short	57445	57449	57447
<b>Long</b>	Small Long	58415	58419	58417
	Medium Long	58425	58429	58427
	Large Long	58435	58439	58437
	X-Large Long	58445	58449	58447



**Ease Opaque Thigh Highs**

Thigh highs in a men's and women's style feature a silicone band for an improved grip and comfortable, all-day wear. A reinforced toe and heel hold up longer.



**15-20mmHg\* Mild Support\*\***

	Women's Black	Women's Sand	Women's Natural	Women's Navy	Men's Black	Men's Khaki	
<b>Short</b>	Small Short	52415	52410	52413	52417	52915	52919
	Medium Short	52425	52420	52423	52427	52925	52929
	Large Short	52435	52430	52433	52437	52935	52939
	X-Large Short					52945	52949
<b>Long</b>	Small Long	53415	53410	53413	53417	53915	53919
	Medium Long	53425	53420	53423	53427	53925	53929
	Large Long	53435	53430	53433	53437	53935	53939
	X-Large Long					53945	53949

**20-30mmHg\* Moderate Support\*\***

	Women's Black	Women's Sand	Women's Natural	Women's Navy	Men's Black	Men's Khaki	
<b>Short</b>	Small Short	50215	50210	50213	50217	51615	51619
	Medium Short	50225	50220	50223	50227	51625	51629
	Large Short	50235	50230	50233	50237	51635	51639
	X-Large Short					51645	51649
<b>Long</b>	Small Long	54215	54210	54213	54217	55615	55619
	Medium Long	54225	54220	54223	54227	55625	55629
	Large Long	54235	54230	54233	54237	55635	55639
	X-Large Long					55645	55649

**30-40mmHg\* Firm Support\*\***

	Women's Black	Women's Sand	Men's Black	Men's Khaki	
<b>Short</b>	Small Short	56215	56210	57615	57619
	Medium Short	56225	56220	57625	57629
	Large Short	56235	56230	57635	57639
	X-Large Short			57645	57649
<b>Long</b>	Small Long	59215	59210	58615	58619
	Medium Long	59225	59220	58625	58629
	Large Long	59235	59230	58635	58639
	X-Large Long			58645	58649



## Ease Opaque Open-Toe Thigh Highs

Open-toe thigh highs combine everyday comfort with a stay-put silicone band for improved grip. An open-toe construction pair well with open-toe footwear, while the reinforced toe and heel offer durability.



Black Sand

### 15–20mmHg\* Mild Support\*\*


	Black	Sand
● Short	Small Short	51315 51310
	Medium Short	51325 51320
	Large Short	51335 51330
	X-Large Short	51345 51340
● Long	Small Long	51815 51810
	Medium Long	51825 51820
	Large Long	51835 51830
	X-Large Long	51845 51840

### 20–30mmHg\* Moderate Support\*\*

	Black	Sand
● Short	Small Short	50715 50710
	Medium Short	50725 50720
	Large Short	50735 50730
	X-Large Short	50745 50740
● Long	Small Long	54715 54710
	Medium Long	54725 54720
	Large Long	54735 54730
	X-Large Long	54745 54740

### 30–40mmHg\* Firm Support\*\*

	Black	Sand
● Short	Small Short	56715 56710
	Medium Short	56725 56720
	Large Short	56735 56730
	X-Large Short	56745 56740
● Long	Small Long	59715 59710
	Medium Long	59725 59720
	Large Long	59735 59730
	X-Large Long	59745 59740



## Ease Opaque Pantyhose

Pantyhose feature a figure flattering control top and a soft, non-restrictive knit-in waistband for comfortable, all-day wear. A reinforced toe and heel hold up to everyday wear and tear.



Black Sand Natural Navy

See size charts on page 14.

**15-20mmHg\* Mild Support\*\***


		Black	Sand	Natural	Navy
Short	Small Short	52615	52610	52613	52617
	Medium Short	52625	52620	52623	52627
	Large Short	52635	52630	52633	52637
	X-Large Short	52645	52640	52643	52647
Long	Small Long	53615	53610	53613	53617
	Medium Long	53625	53620	53623	53627
	Large Long	53635	53630	53633	53637
	X-Large Long	53645	53640	53643	53647

**20-30mmHg\* Moderate Support\*\***

		Black	Sand	Natural	Navy
Short	Small Short	50415	50410	50413	50417
	Medium Short	50425	50420	50423	50427
	Large Short	50435	50430	50433	50437
	X-Large Short	50445	50440	50443	50447
Long	Small Long	54415	54410	54413	54417
	Medium Long	54425	54420	54423	54427
	Large Long	54435	54430	54433	54437
	X-Large Long	54445	54440	54443	54447

**30-40mmHg\* Firm Support\*\***

		Black	Sand
Short	Small Short	56415	56410
	Medium Short	56425	56420
	Large Short	56435	56430
	X-Large Short	56445	56440
Long	Small Long	59415	59410
	Medium Long	59425	59420
	Large Long	59435	59430
	X-Large Long	59445	59440



**Ease Opaque  
Open-Toe Waist Highs**

Open-toe waist highs offer a smooth top and a gentle, knit-in waistband that stays comfortable throughout the day. An open-toe construction pairs easily with open-toe shoes, and a reinforced toe and heel enhance lasting wear.

**MATERNITY-APPROVED**



**15-20mmHg\* Mild Support\*\***

		Black	Sand
Short	Small Short	58715	58710
	Medium Short	58725	58720
	Large Short	58735	58730
Long	Small Long	51915	51910
	Medium Long	51925	51920
	Large Long	51935	51930

**20-30mmHg\* Moderate Support\*\***

		Black	Sand
Short	Small Short	50815	50810
	Medium Short	50825	50820
	Large Short	50835	50830
Long	Small Long	54815	54810
	Medium Long	54825	54820
	Large Long	54835	54830

**30-40mmHg\* Firm Support\*\***

		Black	Sand
Short	Small Short	56815	56810
	Medium Short	56825	56820
	Large Short	56835	56830
Long	Small Long	59815	59810
	Medium Long	59825	59820
	Large Long	59835	59830



## Ease Sheer Knee Highs

These opaque knee highs offer all-day comfort with a soft, knit-in band and a refined sheer-knit texture.

Reinforced toe and heel details ensure lasting durability.



Black Sand Natural Navy

### 15–20mmHg\* Mild Support\*\*

	Black	Sand	Natural	Navy	
Short	Small Short	60015	60010	60013	60017
	Medium Short	60025	60020	60023	60027
	Large Short	60035	60030	60033	60037
	X-Large Short	60045	60040	60043	
Long	Small Long	62415	62410	62413	62417
	Medium Long	62425	62420	62423	62427
	Large Long	62435	62430	62433	62437
	X-Large Long	62445	62440	62443	

### 20–30mmHg\* Moderate Support\*\*

	Black	Sand	Natural	Navy	
Short	Small Short	64015	64010	64013	64017
	Medium Short	64025	64020	64023	64027
	Large Short	64035	64030	64033	64037
	X-Large Short	64045	64040	64043	
Long	Small Long	65415	65410	65413	65417
	Medium Long	65425	65420	65423	65427
	Large Long	65435	65430	65433	65437
	X-Large Long	65445	65440	65443	

### 30–40mmHg\* Firm Support\*\*

	Black	Sand	Natural	
Short	Small Short	66015	66010	66013
	Medium Short	66025	66020	66023
	Large Short	66035	66030	66033
	X-Large Short	66045	66040	66043
Long	Small Long	69415	69410	69413
	Medium Long	69425	69420	69423
	Large Long	69435	69430	69433
	X-Large Long	69445	69440	69443



## Ease Sheer

### Open-Toe Knee Highs

Open-toe knee highs are crafted for everyday wear, featuring a smooth, knit-in band for a comfortable fit that stays put. The open-toe design makes them ideal for open-toe footwear.



Black Sand Natural

See size charts on page 14.

**15-20mmHg\* Mild Support\*\***

		<b>Black</b>	<b>Sand</b>	<b>Natural</b>
<b>Short</b>	Small Short	60115	60110	60113
	Medium Short	60125	60120	60123
	Large Short	60135	60130	60133
<b>Long</b>	Small Long	62515	62510	62513
	Medium Long	62525	62520	62523
	Large Long	62535	62530	62533

**20-30mmHg\* Moderate Support\*\***

		<b>Black</b>	<b>Sand</b>	<b>Natural</b>
<b>Short</b>	Small Short	64115	64110	64113
	Medium Short	64125	64120	64123
	Large Short	64135	64130	64133
<b>Long</b>	Small Long	65515	65510	65513
	Medium Long	65525	65520	65523
	Large Long	65535	65530	65533



**Ease Sheer  
Thigh Highs**

Sheer thigh highs feature a silicone band for an improved grip and comfortable, all-day wear. A reinforced toe and heel hold up longer.



Black Sand Natural Navy

**15-20mmHg\* Mild Support\*\***

		<b>Black</b>	<b>Sand</b>	<b>Natural</b>	<b>Navy</b>
<b>Short</b>	Small Short	60215	60210	60213	60217
	Medium Short	60225	60220	60223	60227
	Large Short	60235	60230	60233	60237
<b>Long</b>	Small Long	62615	62610	62613	62617
	Medium Long	62625	62620	62623	62627
	Large Long	62635	62630	62633	62637

**20-30mmHg\* Moderate Support\*\***

		<b>Black</b>	<b>Sand</b>	<b>Natural</b>	<b>Navy</b>
<b>Short</b>	Small Short	64215	64210	64213	64217
	Medium Short	64225	64220	64223	64227
	Large Short	64235	64230	64233	64237
<b>Long</b>	Small Long	65615	65610	65613	65617
	Medium Long	65625	65620	65623	65627
	Large Long	65635	65630	65633	65637

**30-40mmHg\* Firm Support\*\***

		<b>Black</b>	<b>Sand</b>	<b>Natural</b>
<b>Short</b>	Small Short	66215	66210	66213
	Medium Short	66225	66220	66223
	Large Short	66235	66230	66233
<b>Long</b>	Small Long	69615	69610	69613
	Medium Long	69625	69620	69623
	Large Long	69635	69630	69633



## Ease Sheer

### Open-Toe Thigh Highs

Sheer open-toe thigh highs feature a stay-put silicone band for improved grip. An open-toe construction pair well with open-toe footwear, while the reinforced toe and heel offer durability.



Black Sand Natural

### 15–20mmHg\* Mild Support\*\*

		Black	Sand	Natural
Short	Small Short	60315	60310	60313
	Medium Short	60325	60320	60323
	Large Short	60335	60330	60333
Long	Small Long	62715	62710	62713
	Medium Long	62725	62720	62723
	Large Long	62735	62730	62733

### 20–30mmHg\* Moderate Support\*\*

		Black	Sand	Natural
Short	Small Short	64315	64310	64313
	Medium Short	64325	64320	64323
	Large Short	64335	64330	64333
Long	Small Long	65715	65710	65713
	Medium Long	65725	65720	65723
	Large Long	65735	65730	65733



## Ease Sheer

### Pantyhose

Sheer pantyhose feature a soft knit-in waistband to prevent rolling and enhances comfort. A reinforced toe and heel hold up to everyday wear and tear.



Black Sand Natural Navy

See size charts on page 14.

**15-20mmHg\* Mild Support\*\***

		<b>Black</b>	<b>Sand</b>	<b>Natural</b>	<b>Navy</b>
<b>Short</b>	Small Short	60415	60410	60413	60417
	Medium Short	60425	60420	60423	60427
	Large Short	60435	60430	60433	60437
<b>Long</b>	Small Long	62815	62810	62813	62817
	Medium Long	62825	62820	62823	62827
	Large Long	62835	62830	62833	62837

**20-30mmHg\* Moderate Support\*\***

		<b>Black</b>	<b>Sand</b>	<b>Natural</b>	<b>Navy</b>
<b>Short</b>	Small Short	64415	64410	64413	64417
	Medium Short	64425	64420	64423	64427
	Large Short	64435	64430	64433	64437
<b>Long</b>	Small Long	65815	65810	65813	65817
	Medium Long	65825	65820	65823	65827
	Large Long	65835	65830	65833	65837

**30-40mmHg\* Firm Support\*\***

		<b>Black</b>	<b>Sand</b>	<b>Natural</b>
<b>Short</b>	Small Short	66415	66410	66413
	Medium Short	66425	66420	66423
	Large Short	66435	66430	66433
<b>Long</b>	Small Long	69815	69810	69813
	Medium Long	69825	69820	69823
	Large Long	69835	69830	69833



## Ease Sheer

### Open-Toe Pantyhose

Sheer open-toe pantyhose feature a knit-in waistband that stays comfortable throughout the day. An open-toe construction pairs easily with open-toe shoes, and a reinforced toe and heel enhance lasting wear.



Black Sand Natural

**15-20mmHg\* Mild Support\*\***

		<b>Black</b>	<b>Sand</b>	<b>Natural</b>
<b>Short</b>	Small Short	60515	60510	60513
	Medium Short	60525	60520	60523
	Large Short	60535	60530	60533
<b>Long</b>	Small Long	62915	62910	62913
	Medium Long	62925	62920	62923
	Large Long	62935	62930	62933

**20-30mmHg\* Moderate Support\*\***

		<b>Black</b>	<b>Sand</b>	<b>Natural</b>
<b>Short</b>	Small Short	64515	64510	64513
	Medium Short	64525	64520	64523
	Large Short	64535	64530	64533
<b>Long</b>	Small Long	65915	65910	65913
	Medium Long	65925	65920	65923
	Large Long	65935	65930	65933



## Ease Bold

### Socks

Versatile compression socks that deliver comfort in bold and colorful patterns. A reinforced toe and heel hold up to everyday wear and tear.



Deep Blue   Evergreen   Natural Stripe   Daisy   Starry Night

#### 15–20mmHg\* Mild Support\*\*

	Deep Blue	Evergreen	Natural Stripe	Daisy	Starry Night
Small	55910	55911	55912	55914	55913
Medium	55920	55921	55922	55924	55923
Large	55930	55931	55932	55934	55933
X-Large	55940	55941	55942	55944	55943

#### 20–30mmHg\* Moderate Support\*\*

	Deep Blue	Evergreen	Natural Stripe	Daisy	Starry Night
Small	57010	57011	57012	57014	57013
Medium	57020	57021	57022	57024	57023
Large	57030	57031	57032	57034	57033
X-Large	57040	57041	57042	57044	57043



## Ease Lymphedema

### Arm Sleeves

Arm sleeves feature a specialized knit pattern at the inside of the elbow bend to help the sleeve stay in place. A soft, silicone-dotted band prevents rolling and provides comfort.



Black   Sand

#### 20–30mmHg\* Moderate Support\*\*

	Black	Sand
Regular	Small	63015 63010
	Medium	63025 63020
	Large	63035 63030
	X-Large	63045 63040
Long	Small Long	63215 63210
	Medium Long	63225 63220
	Large Long	63235 63230
	X-Large Long	63245 63240

#### 30–40mmHg\* Firm Support\*\*

	Black	Sand
Regular	Small	63315 63310
	Medium	63325 63320
	Large	63335 63330
	X-Large	63345 63340
Long	Small Long	63415 63410
	Medium Long	63425 63420
	Large Long	63435 63430
	X-Large Long	63445 63440

See size charts on page 14.



## Ease Lymphedema

### Gauntlets

Lymphedema gauntlets feature ultra-stretchy yarns and seamless engineering for all-day wear. A soft, comfortable non-restrictive wrist provides more comfort.



Black Sand

#### 20-30mmHg\* Moderate Support\*\*

	Black	Sand
Small	63715	63710
Medium	63725	63720
Large	63735	63730
X-Large	63745	63740

#### 30-40mmHg\* Firm Support\*\*

	Black	Sand
Small	63815	63810
Medium	63825	63820
Large	63835	63830
X-Large	63845	63840



## Ease Lymphedema

### Gloves

Lymphedema gloves feature seamless compressive fingers and a soft, comfortable non-restrictive wrist. Ultra-stretchy yarns more comfort and ease of donning.



Black Sand

#### 20-30mmHg\* Moderate Support\*\*

	Black	Sand
Small	63915	63910
Medium	63925	63920
Large	63935	63930
X-Large	63945	63940



Scan the QR Code



To view additional Thuasne® Lymphedema products.

Or, ask your sales representative for more information.

core-spun  
BY THERAFIRM®

## GRADIENT COMPRESSION SOCKS

Core-Spun by Therafirm gradient compression socks look and feel like soft, comfortable everyday socks with the added benefit of true gradient compression. Ultra-stretchy yarns make the socks easy to put on and sock thickness makes them more comfortable for all-day wear. Moisture wicking yarns help wick moisture away from the skin to the exterior of the sock making them cooler and more comfortable.

Core-Spun yarns are made by twisting fibers around an inner fiber, creating one yarn. This process is used to enhance functional properties of fabrics, especially stretch and comfort as well as durability.

Patent # 7,895,863

### Compression Levels Available

10–15mmHg\*  
Light Support

Improves  
circulation

15–20mmHg\*  
Mild Support

Prevents  
mild swelling

20–30mmHg\*  
Moderate Support

Prevents mild  
to moderate  
swelling

30–40mmHg\*  
Firm Support

Prevents  
moderate  
swelling

Refer to page 6–7 for more information on compression levels.



Crew socks feature  
a mesh instep and  
gentle arch support



Knee highs are  
available in fashionable  
colors and patterns



Comfort top is non-restrictive while staying in place all day

Ultra-stretchy Core-Spun yarns make sock easy to don

Gentle compression squeeze helps to support muscles and tendons

Sock thickness provides a casual and comfortable alternative to sheer compression

Superior moisture wicking fibers

Delivers lab tested true gradient compression to promote better blood flow, prevent mild to moderate swelling, and relieve tired, achy legs and feet

Smooth toe seam

Reinforced heel and toe



Cushioned Support Socks include a sole with added thickness for extra comfort



Thigh high length available

Ease

Core-Spun

Core-Sport


TheraSport

Therafirm

Compression Solutions


## Core-Spun Size Charts

### Crew Socks



Size	Ankle	Calf	Length
Small	6.5"-8.5" (17-22cm)	11"-16.5" (28-42cm)	12"-16" (30-41cm)
Medium	8"-10" (21-26cm)	12"-17.5" (31-45cm)	13"-17" (33-43cm)
Large	9"-11.5" (23-29cm)	13"-19" (33-49cm)	14"-18" (36-46cm)
X-Large	11"-15" (28-38cm)	17"-23" (43-58cm)	15"-20" (38-51cm)
XX-Large	13"-17" (33-43cm)	18.5"-25.5" (47-65cm)	16"-21" (41-53cm)


### Knee High Socks



Size	Ankle	Calf	Length
Small	6.5"-8.5" (17-22cm)	11"-16.5" (28-42cm)	Short 16" (41cm) or less
Medium	8"-10" (21-26cm)	12"-17.5" (31-45cm)	
Large	9"-11.5" (23-29cm)	13"-19" (33-49cm)	..... Regular Greater than 16" (41 cm))
X-Large	11"-15" (28-38cm)	17"-23" (43-58cm)	
XX-Large	13"-17" (33-43cm)	18.5"-25.5" (47-65cm)	

**Fits calf up to 25.5"**

### Thigh High Socks



Size	Ankle	Calf	Thigh	Length
Small	6.5"-8.5" (17-22cm)	11"-16.5" (28-42cm)	15.5"-21.5" (40-54cm)	Short 28.5" (72.5cm) or less
Medium	8"-10" (21-26cm)	12"-17.5" (31-45cm)	17.5"-25.5" (45-65cm)	
Large	9"-11.5" (23-29cm)	13"-19" (33-49cm)	21.5"-29.5" (55-75cm)	..... Long Greater than 28.5" (72.5cm)
X-Large	11"-15" (28-38cm)	17"-23" (43-58cm)	25.5"-33.5" (65-85cm)	
XX-Large	13"-17" (33-43cm)	18.5"-25.5" (47-65cm)	27.5"-37" (70-94cm)	
XXX-Large	14"-19" (36-48cm)	20.5"-27.5" (52-70cm)	29.5"-39.5" (75-100cm)	

**Fits thigh up to 39.5"**



## Core-Spun

### Crew Socks

Core-Spun crew socks feature an athletic style and light compression. A mesh instep and gentle arch support provide comfort and breathability.



Black White

### 10–15mmHg\* Light Support\*\*

	<b>Black</b>	<b>White</b>
Small	16213	16013
Medium	16223	16023
Large	16233	16033
X-Large	16243	16043
XX-Large	16253	16053



## Core-Spun

### Cushioned Knee High Socks

Engineered with ultra-stretchy Core-Spun yarns for more comfort and ease of donning. Additional thickness in the sole of the sock provide more cushioning and comfort.



Black White

### 15–20mmHg\* Mild Support †

	<b>Black</b>	<b>White</b>
Small	17912	17312
Medium	17922	17322
Large	17932	17332
X-Large	17942	17342
XX-Large	17952	17352

### 20–30mmHg\* Moderate Support\*\*

	<b>Black</b>	<b>White</b>
Small	18912	18312
Medium	18922	18322
Large	18932	18332
X-Large	18942	18342
XX-Large	18952	18352

Ease

Core-Spun

Core-Sport

TheraSport

Therafirm

Compression Solutions




## Core-Spun

### Knee High Socks

Core-Spun knee high socks are knit with ultra-stretchy Core-Spun fibers making to provide more comfort and easy donning. Available in several neutral colors and colorful patterns.



Black White Navy Khaki



Thin Line Merger Trend-setter



Classic Diamond Pink Argyle Slate Argyle Mono-gradient

### 10–15mmHg\* Light Support\*\*

	Black	White	Navy	Khaki	Thin Line	Trend-setter	Slate Argyle
<b>Short</b>	Small Short	16711	16111	16411	16511		
	Medium Short	16721	16121	16421	16521		
	Large Short	16731	16131	16431	16531		
	X-Large Short	16741	16141	16441	16541		
	XX-Large Short	16751	16151	16451	16551		
<b>Regular</b>	Small	16712	16112	16412	16512	11662	11612 12612
	Medium	16722	16122	16422	16522	11672	11622 12622
	Large	16732	16132	16432	16532	11682	11632 12632
	X-Large	16742	16142	16442	16542	11692	11642 12642
	XX-Large	16752	16152	16452	16552		

### 30–40mmHg\* Firm Support\*\*

	Black	White	Khaki
<b>Short</b>	Small Short	19711	19111 19511
	Medium Short	19721	19121 19521
	Large Short	19731	19131 19531
	X-Large Short	19741	19141 19541
	XX-Large Short	19751	19151 19551
<b>Regular</b>	Small	19712	19112 19512
	Medium	19722	19122 19522
	Large	19732	19132 19532
	X-Large	19742	19142 19542
	XX-Large	19752	19152 19552

### 15–20mmHg\* Mild Support†

	Black	White	Navy	Khaki	Thin Line	Merger	Trend-setter	Classic Diamond	Pink Argyle	Slate Argyle	Mono-gradient
<b>Short</b>	Small Short	17711	17111	17411	17511						
	Medium Short	17721	17121	17421	17521						
	Large Short	17731	17131	17431	17531						
	X-Large Short	17741	17141	17441	17541						
	XX-Large Short	17751	17151	17451	17551						
<b>Regular</b>	Small	17712	17112	17412	17512	11762	10762	11712	10712	12762	12712 13712
	Medium	17722	17122	17422	17522	11772	10772	11722	10722	12772	12722 13722
	Large	17732	17132	17432	17532	11782	10782	11732	10732	12782	12732 13732
	X-Large	17742	17142	17442	17542	11792	10792	11742	10742	12792	12742 13742
	XX-Large	17752	17152	17452	17552						

### 20–30mmHg\* Moderate Support\*\*

	Black	White	Navy	Khaki	Thin Line	Merger	Trend-setter	Classic Diamond	Pink Argyle	Slate Argyle	Mono-gradient
<b>Short</b>	Small Short	18711	18111	18411	18511						
	Medium Short	18721	18121	18421	18521						
	Large Short	18731	18131	18431	18531						
	X-Large Short	18741	18141	18441	18541						
	XX-Large Short	18751	18151	18451	18551						
<b>Regular</b>	Small	18712	18112	18412	18512	11862	10862	11812	10812	12862	12812 13812
	Medium	18722	18122	18422	18522	11872	10872	11822	10822	12872	12822 13822
	Large	18732	18132	18432	18532	11882	10882	11832	10832	12882	12832 13832
	X-Large	18742	18142	18442	18542	11892	10892	11842	10842	12892	12842 13842
	XX-Large	18752	18152	18452	18552						


See size charts on page 28.



## Core-Spun Thigh High Socks

Full-length compression socks engineered with ultra-stretchy Core-Spun fibers for more comfort and easy donning. A comfortable, non-restrictive silicone band that stays up all day.

**FITS THIGH UP TO 39.5"**



Black White

### 15-20mmHg\* Mild Support†

	Black	White
● Short	Small Short	17617 17117
	Medium Short	17627 17127
	Large Short	17637 17137
	X-Large Short	17647 17147
	XX-Large Short	17657 17157
	XXX-Large Short	17667 17167
● Long	Small Long	17618 17118
	Medium Long	17628 17128
	Large Long	17638 17138
	X-Large Long	17648 17148
	XX-Large Long	17658 17158
	XXX-Large Long	17668 17168

### 20-30mmHg\* Moderate Support\*\*

	Black	White
● Short	Small Short	18617 18117
	Medium Short	18627 18127
	Large Short	18637 18137
	X-Large Short	18647 18147
	XX-Large Short	18657 18157
	XXX-Large Short	18667 18167
● Long	Small Long	18618 18118
	Medium Long	18628 18128
	Large Long	18638 18138
	X-Large Long	18648 18148
	XX-Large Long	18658 18158
	XXX-Large Long	18668 18168

### 30-40mmHg\* Firm Support\*\*

	Black	
● Short	Small Short	19617
	Medium Short	19627
	Large Short	19637
● Long	Small Long	19618
	Medium Long	19628
	Large Long	19638



## ATHLETIC GRADIENT COMPRESSION SOCKS

Core-Sport by Therafirm athletic compression socks provide a controlled amount of pressure which is greatest at the ankle and gradually decreases toward the top. True gradient compression helps to improve circulation, which can increase energy for endurance, improve performance, reduce muscle fatigue, and speed up recovery time after a workout. Core-Sport socks feature soft, ultra-stretchy Core-Spun yarns for a comfortable performance sock.

Patent #7,895,863

### Compression Levels Available

**15–20mmHg\***  
Mild Support


**20–30mmHg\***  
Moderate Support

Promotes faster recovery

Enhances performance

Refer to page 6–7 for more information on compression levels.

### Size Chart



Size	Ankle	Calf	Length
Small	6.5"–8.5" (17–22cm)	11"–16.5" (28–42cm)	12"–16" (30–41cm)
Medium	8"–10" (21–26cm)	12"–17.5" (31–45cm)	13"–17" (33–43cm)
Large	9"–11.5" (23–29cm)	13"–19" (33–49cm)	14–18" (36–46cm)
X-Large	11"–15" (28–38cm)	17"–23" (43–58cm)	15–20" (38–51cm)



Comfort top stays in place without binding

Featuring patented Core-Spun yarn technology

Graduated compression for better circulation and added support and stabilization

Advanced, technical fibers wick moisture away from the skin to the exterior of the sock

Arch compression for added support

Antimicrobial to help prevent odor in the sock

Cushioned heel and foot reduce impact



Black    Grey    White

## Core-Sport Athletic Socks

**15-20mmHg\*** Mild Support\*\*

	<b>Black</b>	<b>Grey</b>	<b>White</b>
Small	37712	37612	37512
Medium	37722	37622	37522
Large	37732	37632	37532
X-Large	37742	37642	37542

**20-30mmHg\*** Moderate Support\*\*

	<b>Black</b>
Small	87712
Medium	87722
Large	87732
X-Large	87742

Ease

Core-Spun

Core-Sport

TheraSport

Therafirm

Compression Solutions

**THERASPORT**  
BY THERAFIRM®

**ATHLETIC GRADIENT  
COMPRESSION SOCKS +  
LEG SLEEVES**

TheraSport by Therafirm athletic compression socks and leg sleeves feature true gradient compression, which delivers a controlled amount of pressure that is greatest at the ankle and gradually decreases as it comes up the leg. Gradient compression helps improve circulation, which can help energize legs, help reduce fatigue, and speed post-workout recovery. TheraSport socks feature comfortable, lightweight fibers for a cooler feel.

**Compression Levels Available**

**15–20mmHg\***  
Mild Support


**20–30mmHg\***  
Moderate Support

Promotes faster recovery

Enhances performance

Refer to page 6–7 for more information on compression levels.

**Size Chart**



Size	Ankle	Calf	Length
Small	6.5"–8.5" (17–22cm)	11"–16.5" (28–42cm)	12"–16" (30–41cm)
Medium	8"–10" (21–26cm)	12"–17.5" (31–45cm)	13"–17" (33–43cm)
Large	9"–11.5" (23–29cm)	13"–19" (33–49cm)	14–18" (36–46cm)
X-Large	11"–15" (28–38cm)	17"–23" (43–58cm)	15–20" (38–51cm)



Non-binding 1" comfort band that stays up

High tech moisture wicking yarns

Mesh ventilation for breathability

"Y" stitch heel that prevents slide



Wear any socks with Leg Sleeves that feature footless design

Cushioning absorbs shock and protects against it

Arch band provides metatarsal support



Black White Pink Blue Red

## TheraSport

### Athletic Socks

**15-20mmHg\*** Mild Support\*\*

	Black	White	Pink	Blue	Red
Small	37415	37413	37410	37417	37416
Medium	37425	37423	37420	37427	37426
Large	37435	37433	37430	37437	37436
X-Large	37445	37443	37440	37447	37446

**20-30mmHg\*** Moderate Support\*\*

	Black	White	Pink	Blue	Red
Small	87415	87413	87410	87417	87416
Medium	87425	87423	87420	87427	87426
Large	87435	87433	87430	87437	87436
X-Large	87445	87443	87440	87447	87446

## TheraSport

### Athletic Leg Sleeves

**15-20mmHg\*** Mild Support\*\*

	Black	White	Pink	Blue
Small	67415	67413	67410	67417
Medium	67425	67423	67420	67427
Large	67435	67433	67430	67437
X-Large	67445	67443	67440	67447

**20-30mmHg\*** Moderate Support\*\*

	Black	White	Pink	Blue
Small	77415	77413	77410	77417
Medium	77425	77423	77420	77427
Large	77435	77433	77430	77437
X-Large	77445	77443	77440	77447

Ease

Core-Spun

Core-Sport

TheraSport

Therafirm

Compression Solutions

# THERAFIRM®

## GRADIENT COMPRESSION HOSIERY

Therafirm compression products are soft, cool and comfortable for all-day wear. Therafirm gradient compression hosiery garments are made using moisture wicking fibers, which cause moisture to evaporate quickly and provide a comfortable coolness. Engineered with advanced technologies to deliver health benefits.

### Compression Levels Available

**10–15mmHg\***  
Light Support

Improves  
circulation

**15–20mmHg\***  
Mild Support

Prevents mild  
swelling


**20–30mmHg\***  
Moderate Support

Prevents mild  
to moderate  
swelling

**30–40mmHg\***  
Firm Support

Prevents  
moderate  
swelling

Refer to page 6–7 for more information on compression levels.



Trouser socks feature  
a wide ribbed style  
and everyday colors




Thigh highs feature a  
silicone dot band that keeps  
them in place all day



Pantyhose feature a non-restrictive waistband that prevents rolling and helps ensure a comfortable fit for all day wear



Open-toe styles available



Soft comfort band is non-restrictive while staying in place all day

Fashionable semi-opaque knit

Two-way stretch nylon and spandex blend

Moisture wicking fibers provide a comfortable coolness

Lab tested true gradient compression to promote better blood flow, prevent mild to moderate swelling, and relieve tired, achy legs and feet

Reinforced toe for durability

Ease

Core-Spun

Core-Sport

TheraSport


Therafirm

Compression Solutions


## Therafirm Size Charts

### Knee Highs

(10–15mmHg\*)



Size	Ankle	Calf	Length	Shoe
Small	7"–8" (18–20cm)	11"–15" (28–38cm)	14"–16" (36–41cm)	5–7
Medium	8"–9.5" (20–24cm)	13"–16" (33–41cm)	15"–17" (38–43cm)	7–9
Large	9.5"–11" (24–28cm)	15"–18" (38–46cm)	16"–18" (41–46cm)	9–10.5
X-Large	11"–12.5" (28–32cm)	17.5"–20.5" (42–52cm)	17.5"–19.5" (45–50cm)	10.5–12



### Women's Trouser Socks


(10–15mmHg\*)

Size	Ankle	Calf	Length	Shoe
Small	7.5"–9" (19–22cm)	13"–15.5" (33–39cm)	13"–19" (33–48cm)	5.5–7.5
Medium	9.25"–10.75" (24–27cm)	15"–17.5" (38–44cm)	14"–20" (36–51cm)	7–10
Large	9.5"–11" (24–28cm)	17.75"–18.5" (45–47cm)	14"–18" (36–46cm)	8–11

### Men's Trouser Socks


(10–15mmHg\*)

Size	Ankle	Calf	Length	Shoe
Small	9"–11" (23–28cm)	11"–16.5" (28–42cm)	12"–15" (31–38cm)	7.5–10
Medium	10"–14" (25–36cm)	12"–17" (31–43cm)	13"–16" (33–41cm)	10.5–12
Large	11"–15" (28–38cm)	13"–19" (33–48cm)	14"–18" (36–46cm)	12.5+




### Knee Highs

(15–20mmHg\*)



Size	Ankle	Calf	Length
Small	6.25"–7.25" (16–18cm)	11"–13.5" (28–34cm)	12"–18" (31–46cm)
Medium	7.5"–9" (19–22cm)	13"–15.5" (33–39cm)	13"–19" (33–48cm)
Large	9.25"–10.75" (24–27cm)	15"–17.5" (38–44cm)	14"–20" (36–51cm)
X-Large	11"–11.5" (28–29cm)	17.75"–18.5" (45–47cm)	15"–22" (38–56cm)
XX-Large	11.75"–12.75" (30–32cm)	18"–19.5" (46–50cm)	15"–22" (38–56cm)



### Women's Trouser Socks


(15–20mmHg\* + 20–30mmHg\*)

Size	Ankle	Calf	Length
Small	6.25"–7.25" (16–18cm)	11"–13.5" (28–34cm)	12"–18" (31–46cm)
Medium	7.5"–9" (19–22cm)	13"–15.5" (33–39cm)	13"–19" (33–48cm)
Large	9.25"–10.75" (24–27cm)	15"–17.5" (38–44cm)	14"–20" (36–51cm)
X-Large	11"–11.5" (28–29cm)	17.75"–18.5" (45–47cm)	15"–22" (38–56cm)

### Men's Trouser Socks


(15–20mmHg\*, 20–30mmHg\*, + 30–40mmHg\*)

Size	Ankle	Calf	Length
Small	7"–9" (18–23cm)	10.5"–15" (27–38cm)	11"–14" (28–36cm)
Medium	9"–11" (23–28cm)	11"–16.5" (28–42cm)	12"–15" (31–38cm)
Large	10"–14" (25–36cm)	12"–17" (31–43cm)	13"–16" (33–41cm)
X-Large	11"–15" (28–38cm)	13"–19" (33–48cm)	14"–18" (36–46cm)



### Knee Highs + Open-Toe Knee Highs

(20–30mmHg\* + 30–40mmHg\*)



Size	Ankle	Calf	Length
Small	6.25"–7.25" (16–18cm)	10.5"–13" (27–33cm)	11"–15" (28–38cm)
Medium	7.5"–9" (19–22cm)	11.5"–14" (29–36cm)	12"–17" (31–43cm)
Large	9"–10" (22–25cm)	12.5"–15" (32–38cm)	13"–18" (33–46cm)
X-Large	10"–11.25" (25–29cm)	13.5"–16" (34–41cm)	14"–19" (36–48cm)
XX-Large	11.25"–12.5" (29–32cm)	14.5"–17" (37–43cm)	16"–21" (41–53cm)
XXX-Large	12.5"–13.75" (32–35cm)	16.5"–20.5" (42–52cm)	16"–21" (41–53cm)
XXXX-Large	13.75"–15" (35–38cm)	20"–24" (51–61cm)	16"–21" (41–53cm)

## Thigh Highs


(15–20mmHg\*, 20–30mmHg\*, + 30–40mmHg\*)

## Thigh Highs with Lace Top

(20–30mmHg\*)

## Open-Toe Thigh Highs


(15–20mmHg\* + 20–30mmHg\*)



Size	Ankle	Calf	Thigh	Length
Small	6.25"–7.5" (16–19cm)	10.5"–13" (27–33cm)	15"–18" (38–46cm)	23"–26" (58–66cm)
Medium	7.5"–8.75" (19–22cm)	11.5"–14" (29–36cm)	17.5"–20.5" (44–52cm)	23"–29" (58–74cm)
Large	8.75"–10" (22–25cm)	12.5"–15" (32–38cm)	20"–23" (51–58cm)	25"–31" (64–79cm)
X-Large	10"–11.25" (25–29cm)	13.5"–16" (34–41cm)	22.5"–25.5" (57–65cm)	27"–33" (69–84cm)
XX-Large	11.25"–12.5" (29–32cm)	14.5"–17" (37–43cm)	25"–28" (64–71cm)	29"–35" (74–89cm)

## Microfiber Footless Tights

Current Weight




## Pantyhose

(20–30mmHg\*)

Size	Ankle	Calf	Thigh	Hip	Waist	Length
Small	6.25"–7.5" (16–19cm)	10.5"–13" (27–33cm)	15"–18" (38–46cm)	to 40" (102cm)	to 30" (76cm)	29"–35" (74–89cm)
Medium	7.5"–8.75" (19–22cm)	11.5"–14" (30–36cm)	17.5"–20.5" (44–52cm)	to 42" (107cm)	to 33.5" (85cm)	31"–37" (79–94cm)
Large	8.75"–10" (22–25cm)	12.5"–15" (32–38cm)	20"–23" (51–58cm)	to 44" (112cm)	to 36" (91cm)	33"–39" (84–99cm)
X-Large	10"–11.25" (25–29cm)	13.5"–16" (34–41cm)	22.5"–25.5" (57–65cm)	to 46" (117cm)	to 38" (97cm)	35"–41" (89–104cm)
XX-Large	11.25"–12.5" (29–32cm)	14.5"–17" (37–43cm)	25"–28" (64–71cm)	to 52" (132cm)	to 44" (112cm)	35"–41" (89–104cm)

## Thigh Highs with Lace Top

(15–20mmHg\*)




Size	Ankle	Calf	Thigh	Length
Small	6.25"–7.25" (16–18cm)	11"–13.5" (28–34cm)	14"–15.5" (36–39cm)	23"–25" (58–64cm)
Medium	7.5"–9" (19–23cm)	13"–15.5" (33–39cm)	16"–17.5" (41–44cm)	25"–28" (64–71cm)
Large	9.25"–10.75" (23–27cm)	15"–17.5" (38–44cm)	17.5"–19" (44–48cm)	28"–31" (71–79cm)
X-Large	11"–11.5" (28–29cm)	17.75"–18.5" (45–47cm)	19"–21" (48–53cm)	31"–34" (79–86cm)
XX-Large	11.75"–12.5" (30–32cm)	18"–19.5" (46–50cm)	21"–23" (53–58cm)	34"–36" (86–91cm)

## Pantyhose

(15–20mmHg\*)

Current Weight



Ease

Core-Spun

Core-Sport

TheraSport

Therafirm

Compression Solutions



## Therafirm Knee Highs

A soft comfort band is non-restrictive while staying in place all day. A reinforced toe and heel provide more durability.



Black Sand Natural White

### 10–15mmHg\* Light Support\*\*

	Black	Sand	Natural	White
Small	33015	33016	33014	33017
Medium	33025	33026	33024	33027
Large	33035	33036	33034	33037
X-Large	33045	33046	33044	33047

### 15–20mmHg\* Mild Support\*\*

	Black	Sand	Natural
Small	68128	68200	68193
Medium	68124	68201	68194
Large	68125	68202	68195
X-Large	68132	68203	68196
XX-Large	68121	68204	68197

### 20–30mmHg\* Moderate Support\*\*

	Black	Sand
Small	77108	75112
Medium	77104	75113
Large	77105	75114
X-Large	77106	75115
XX-Large	77107	75116
XXX-Large	77109	75117
XXXX-Large	77111	75118


### 30–40mmHg\* Firm Support\*\*

	Black	Sand
Small	77600	76605
Medium	77601	76606
Large	77602	76607
X-Large	77603	76608
XX-Large	77604	76609



## Therafirm Open-Toe Knee Highs

Features open-toe style and soft comfort band is non-restrictive while staying in place all day. Two-way stretch nylon and spandex.



Black Sand

### 20–30mmHg\* Moderate Support\*\*

	Black	Sand
Small	77300	79320
Medium	77301	79321
Large	77302	79322
X-Large	77303	79323
XX-Large	77304	79324
XXX-Large	77305	79325
XXXX-Large	77306	79326

### 30–40mmHg\* Firm Support\*\*

	Black	Sand
Small	77500	76505
Medium	77501	76506
Large	77502	76507
X-Large	77503	76508
XX-Large	77504	76509


See size charts on page 38.



## Therafirm

### Trouser Socks

Ribbed trouser socks feature over-the-calf compression with a non-binding, knit-in band to keep socks in place all day. Heel and toe are engineered to provide more comfort.



Black Navy Khaki White Brown

#### 10–15mmHg\* Light Support\*\*

	Women's Black	Men's Black
Small	90215	90415
Medium	90225	90425
Large	90235	90435

#### 15–20mmHg\* Mild Support†

	Women's Black	Women's Khaki	Women's White	Men's Black	Men's Brown	Men's Navy	Men's Khaki	Men's White
Small	68500	68525	68515	68310	68353	68350	68363	68360
Medium	68501	68526	68516	68311	68313	68319	68326	68324
Large	68502	68527	68517	68312	68314	68320	68327	68325
X-Large	68503	68528	68518	68340	68342	68341	68362	68361

#### 20–30mmHg\* Moderate Support\*\*

	Women's Black	Women's Khaki	Men's Black	Men's Brown	Men's Navy	Men's Khaki	Men's White
Small	69500	69525	68315	68347	68346	68349	68331
Medium	69501	69526	68316	68301	68303	68321	68307
Large	69502	69527	68317	68302	68304	68322	68308
X-Large	69503	69528	68318	68328	68309	68323	68330

#### 30–40mmHg\* Firm Support\*\*

	Men's Black
Small	69315
Medium	69316
Large	69317
X-Large	69318



## Therafirm

### Thigh Highs with Lace-Top

Sheer thigh high stockings knit with two-way stretch nylon and spandex feature non-restrictive lace-top band with silicone border to keep stockings in place. A reinforced toe and heel provide more durability.



Black Sand

#### 15–20mmHg\* Mild Support\*\*

	Black	Sand
Small	68415	68450
Medium	68416	68451
Large	68417	68452
X-Large	68418	68453
XX-Large	68419	68454

#### 20–30mmHg\* Moderate Support\*\*

	Black	Sand
Small	71100	71305
Medium	71101	71306
Large	71102	71307
X-Large	71103	71308
XX-Large	71104	71309

Ease


Core-Spun

Core-Sport

TheraSport

Therafirm

Compression Solutions



## Therafirm Thigh Highs

Thigh highs are knit with two-way stretch nylon and spandex and feature a soft, non-restrictive silicone dotted band to keep stockings from slipping. A durable toe and heel are reinforced to last longer.



Black



Sand

### 15–20mmHg\* Mild Support\*\*

	Black	Sand
Small	68465	68460
Medium	68466	68461
Large	68467	68462
X-Large	68468	68463
XX-Large	68469	68464

### 20–30mmHg\* Moderate Support\*\*

	Black	Sand
Small	77200	74212
Medium	77201	74213
Large	77202	74214
X-Large	77203	74215
XX-Large	77204	74216

### 30–40mmHg\* Firm Support\*\*

	Black	Sand
Small	77700	76705
Medium	77701	76706
Large	77702	76707
X-Large	77703	76708
XX-Large	77704	76709



## Therafirm Open-Toe Thigh Highs

Thigh high style stockings feature an open-toe style and a soft, non-restrictive silicone dotted band that stays in place. Knit with two-way stretch nylon and spandex.



Black



Sand


### 20–30mmHg\* Moderate Support\*\*

	Black	Sand
Small	77100	74112
Medium	77101	74113
Large	77102	74114
X-Large	77103	74115
XX-Large	77110	74116

### 30–40mmHg\* Firm Support\*\*

	Black	Sand
Small	77800	76805
Medium	77801	76806
Large	77802	76807
X-Large	77803	76808
XX-Large	77804	76809

See size charts on page 39.



## Therafirm

### Microfiber Footless Tights

Microfiber tights feature a high waist design to provide shaping support and superior comfort. A footless design allows tights to be worn with more footwear options.



Black

#### 10-15mmHg\* Light Support\*\*

	<b>Black</b>
Small	37115
Medium	37125
Large	37135



## Therafirm

### Pantyhose

Knit with sheer, two-way stretch fibers, pantyhose feature a control top and a soft non-restrictive knit-in waistband for comfortable all-day wear. A reinforced toe and heel provide more durability.



Black



Sand



Natural

#### 15-20mmHg\* Mild Support\*\*

	<b>Black</b>	<b>Sand</b>	<b>Natural</b>
Small	68026	68080	68070
Medium	68027	68081	68071
Large	68028	68082	68072
X-Large	68029	68083	68073
XX-Large	68037	68084	68074

#### 20-30mmHg\* Moderate Support\*\*

	<b>Black</b>	<b>Sand</b>
Small	77000	75405
Medium	77001	75406
Large	77002	75407
X-Large	77003	75408
XX-Large	77004	75409

Ease

Core-Spun

Core-Sport

TheraSport

Therafirm

Compression Solutions

# COMPRESSION SOLUTIONS



Black

30–50mmHg\* Firm Support^

		Black
Regular	Small Regular	68615
	Medium Regular	68625
	Large Regular	68635
	X-Large Regular	68645
Tall	Small Tall	68715
	Medium Tall	68725
	Large Tall	68735
	X-Large Tall	68745

## Therafirm

### Ease Adjust Inelastic Compression Below Knee Wrap

Adjustable gradient inelastic compression wrap for venous and lymphatic disease. Designed to provide 30–50mmHg pressure based on wearer needs.

▶ **Liner Socks Included**

A pair of 10–15mmHg Core-Spun by Therafirm® Socks provide foot and ankle compression

▶ **Product Guarantee**

Therafirm provides a 6-month warranty for materials and workmanship

▶ **Lymphedema Management**

Inelastic compression is also applicable for patients with lymphedema

▶ **Day or Night Use**

Inelastic materials offer a high working pressure during the day and a low resting pressure with strong containment at night

▶ **Maximum Durability**

Made with premium grade materials


▶ **Anatomical Design**

Contours to the calf for a better fit, plus seamless straps reduce pressure points

▶ **Easy Application**

Inner sleeve supports the wrap while donning and adjusting straps

### Size Chart



Size	Ankle	Calf	Under Knee	Length
Small	6.5"–9.5" (16–24cm)	9.5"–14" (24–35cm)	10"–14" (25–36cm)	Regular Less than 14" (36cm)
Medium	8"–11" (20–28cm)	12.5"–17" (32–43cm)	12.5"–17" (32–43cm)	
Large	9.5"–12.5" (24–32cm)	16"–20" (40–51cm)	15.5"–19.5" (39–50cm)	Tall More than 14" (36cm)
X-Large	11"–14" (28–36cm)	19"–23" (48–59cm)	18"–22.5" (46–57cm)	



## Therafirm

### Anti-Embolism Stockings

True gradient compression knee high or thigh high stockings feature an inspection toe with compressive support of 18mmHg\* for the pre- or post-operative patient and may be prescribed for the mobility-impaired or non-ambulatory patient. Designed to help prevent edema, leg discomfort, and deep vein thrombosis in individuals subjected to immobility.




White

18mmHg\*

Mild Support †

	White Knee High	White Thigh High
Small	87101	97101
Medium	87102	97102
Large	87103	97103
X-Large	87104	97104
XX-Large	87105	97105

### Size Chart



Size	Knee High Length	Thigh High Length	Ankle	Calf	Thigh
Small	16" or less (41cm or less)	28" or less (71cm or less)	6.5"-8" (17-20cm)	10"-12" (25-31cm)	15"-18" (38-46cm)
Medium	16"-18" (41-46cm)	28"-30" (71-76cm)	7.5"-9" (19-23cm)	11"-13" (28-33cm)	16"-19" (41-48cm)
Large	17"-19" (43-48cm)	30"-32" (74-81cm)	8.5"-10" (22-25cm)	12"-15" (31-38cm)	17"-21" (43-53cm)
X-Large	18"-22" (46-56cm)	30"-34" (76-86cm)	9"-11.50" (23-29cm)	13"-17" (33-43cm)	18"-23" (46-58cm)
XX-Large	17.5"-19.5" (44-50cm)	29"-33" (74-84cm)	9.5"-13" (24-33cm)	14"-18" (36-46cm)	19"-26" (48-66cm)

Ease

Core-Spun

Core-Sport

TheraSport

Therafirm

Compression Solutions



## Knee Interfaces

A knee brace interface with soft, cool, 4-way stretch, or worn alone for compressive support of 30–40mmHg\* at 50% stretch. Moisture-wicking fibers so skin stays cool, dry and comfortable. Supportive compression and conformity gained from elastic and spandex. Multi-directional stretch moves with body during activity. Inhibits odor in the garment. Made with latex-free materials.

*Sizing: Measure at the largest part of the leg to be covered.*



White Black Grey

### 30–40mmHg\* Firm Support\*\*

#### Straight Sleeves

Black (Polyester)	Flat Width Size	Thigh Circumf.	Length
2BCX-04-18	4" (10cm)	11"–14" (28–36cm)	18" (45cm)
2BCX-05-18	5" (13cm)	14"–17" (36–43cm)	18" (45cm)
2BCX-06-18	6" (15cm)	17"–20" (43–50cm)	18" (45cm)
2BCX-07-18	7" (18cm)	20"–23" (50–58cm)	18" (45cm)
2BCX-07-20	7" (18cm)	20"–23" (50–58cm)	20" (50cm)
2BCX-08-20	8" (20cm)	23"–29" (58–74cm)	20" (50cm)
2BCX-09-20	9" (23cm)	26"–34" (66–86cm)	20" (50cm)

### 30–40mmHg\* Firm Support\*\*

#### Tapered Sleeves

White (Anti-Microbial)	Black (Polyester)	Grey (Polyester)	Flat Width Size	Thigh Circumf.	Length
	2BCX-T5-18	2PCX-T5-18	5" x 4"	14"–17" (36–43cm)	18" (45cm)
	2PCS-T6-18	2PCX-T6-18	6" x 4"	17"–20" (43–50cm)	18" (45cm)
	2PCS-T7-18	2PCX-T7-18	7" x 5"	20"–23" (50–58cm)	18" (45cm)
	2PCS-T8-18	2PCX-T8-18	8" x 6"	23"–29" (58–74cm)	18" (45cm)
	2PCS-T9-18	2PCX-T9-18	9" x 7"	26"–34" (66–86cm)	18" (45cm)



Black

## Compressogrip® Sleeves

A unique elastic tubular bandage for wrists, elbows, arms, ankles, knees, legs, and torso that applies compression for comfort and support under braces. Single layer provides medium compression and double layer provides heavy compression at 50% stretch. At 100% stretch the compression is almost doubled. Doubled layers will almost double compression at both 50% and 100% stretch. Applications of more than two layers should be monitored as compression may exceed the limits of desired therapeutic pressure.

**LATEX-FREE**

Size	Circumf.	18" Length	24" Length
3	9"–12" (22.5–30cm)	2ECA-03-18	
4	12"–16" (30–40cm)	2ECA-04-18	2ECA-04-24
5	15"–20" (38–50cm)	2ECA-05-18	2ECA-05-24
6	18"–24" (45–60cm)	2ECA-06-18	2ECA-06-24
8	24"–32" (60–80cm)	2ECA-08-18	2ECA-08-24



## Compressogrip® Tubular Elastic Bandages

Single layer application provides medium compression, at 50% stretch. Double layer will provide approximately double pressure.

**LATEX-FREE**

Do not use more than two layers unless medically recommended.



White Beige Black

Re-order codes size	Fits Circumf. inches/cm	White (Polyester)	Beige (Polyester)	Black (Polyester)
2" (A)	5-7" (12.5-17.5cm)	6S-7628-02	6S-7631-02	
2.5" (B)	7.5-8" (19-20cm)	6S-7628-S2	6S-7631-S2	
2.75" (C)	8.5-8.75" (21-22cm)	6S-7628-T2	6S-7631-T2	
3" (D)	9-10" (23-25.5cm)	6S-7628-03	6S-7631-03	6S-7633-03
3.5" (E)	10.5-11.5" (26.5-29cm)	6S-7628-S3	6S-7631-S3	
4" (F)	12-14" (30.5-35.5cm)	6S-7628-04	6S-7631-04	6S-7633-04
5" (G)	15-20" (38-51cm)	6S-7628-05	6S-7631-05	6S-7633-05
6" (J)	18-24" (46-61cm)	6S-7628-06	6S-7631-06	6S-7633-06
8" (K)	24-32" (61-81cm)	6S-7628-08	6S-7631-08	6S-7633-08
10" (L)	30-40" (76-102cm)	6S-7628-10	6S-7631-10	
12" (M)	36-48" (91-122cm)	6S-6628-12	6S-6631-12	
14" (N)	43-57" (110-145cm)	6S-6628-14		
16" (P)	48-64" (122-163cm)	6S-6628-16		

**ONLINE ORDERING SYSTEM**  
WHOLESALE.THERAFIRM.COM

**Browse the site and add products to your order within each product page or, if you prefer to place an order using only SKUs, use the Quick Order Pad.**



[wholesale.therafirm.com](http://wholesale.therafirm.com)



## TERMS AND CONDITIONS

### Domestic Accounts

#### TERMS:

Net 30 from date of invoice.

Forms of payment include check, money order, cashier's check, wire transfer or upon approval, credit card (Visa, MasterCard, and Discover) subject to a service fee as allowed by law.

Orders will be shipped C.O.D or CREDIT CARD until credit terms are approved upon submission of a credit application and tax exemption certificate. Upon credit approval, open terms will be established. Credit approval remains subject to continued review even after credit is approved. We will not DROP SHIP C.O.D. ORDERS. If C.O.D. shipments are refused, all costs incurred plus a \$10.00 charge will be billed.

Statements of Account are sent the middle of each month.

Our billing month covers invoices, credit memos and payments for such month. Service charges are based on the amount past due and subject to an interest charge of 1.5% per month, effective annual rate of 18% (or the highest rate permitted by applicable law).

Collection Fee: All collection charges and fees incurred in collecting past due accounts will be billed and paid by customer and must be paid in full before shipment of orders are resumed.

Returned Checks: There will be a \$25.00 service charge on all checks returned unpaid.

#### DELIVERY CLAIMS:

No claim for shortages will be allowed unless made within 10 working days of receipt of order.

#### SALES ONLY IN THE US:

Customers shall not directly or indirectly (including, without limitation, via the Internet) market, sell or export the products to customers located outside of the United States or sell, service or promote the products from a location outside of the United States. If customer is a distributor, customers of the distributor may only sell the products to end users in the United States.

#### ONLINE RESELLERS:

Thuasne customers may not re-sell or distribute to anyone that may re-sell any Thuasne products on any online marketplaces, including but not limited to Amazon.com, Wal-Mart.com, eBay.com, without written consent from Thuasne. Current online retailers are not automatically granted permission to sell on online marketplaces. All customers must have written consent from Thuasne to sell or distribute to customers that sell through these specific channels. Thuasne reserves the right to not sell to any customer who violates this policy.

#### RETURNED GOODS POLICY:

Acceptance of returned merchandise is at the discretion of Thuasne. Returned merchandise will not be accepted for credit without a valid return authorization number and the lot number for the product. There is a maximum 45-day time limit in which to return all merchandise, should Thuasne approve the return.

**Note:** The RA# authorizes your company to return the items, however, based on evaluation, credit to be determined.

#### THE FOLLOWING WILL BE REQUIRED FOR ANY RETURN:

\*Full description of the merchandise.

\*Complete reason for the return

The following applies to all returned merchandise that is not deemed defective or that is not the result of a shipping error:

Please DO NOT return any merchandise that has been worn or tried on. Credit or exchange will be issued at the discretion of Thuasne.

\*0% Restock fee if returned within 45 days and no repackaging required.

\*15% Restock fee if repackaging is required.

\*15% Restock fee if product is returned without an RA.

\*15% Restock fee if returned over 45-days with approval.

#### PAYMENT/ REMITTANCE

120 Osage Avenue  
Kansas City, KS 66105

#### RETURNS

167 Marks Creek Lane  
Hamlet, NC 28345

#### ALL RETURNS MUST BE SENT TO:

167 Marks Creek Lane Hamlet, NC 28345

#### EXPORT ACCOUNTS

See your Sales Representative for specific Terms and Conditions.

#### THUASNE MANUFACTURER'S WARRANTY:

All Thuasne products are warranted to be free of material or workmanship defects for 45 days from date of purchase. At the discretion of Thuasne, any defective merchandise will be replaced or credited at the actual price paid. No refunds.

THIS WARRANTY IS MADE EXPRESSLY IN LIEU OF ANY OTHER WARRANTY AND THUASNE MAKES NO OTHER REPRESENTATION OR WARRANTY OF ANY OTHER KIND, EXPRESS OR IMPLIED, WITH RESPECT TO THE PRODUCTS, WHETHER AS TO MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE INCLUDING ANY PARTICULAR MEDICAL CONDITION, NON-INFRINGEMENT OF THIRD PARTY RIGHTS OR ANY OTHER MATTER. THE WARRANTIES CONTAINED HEREIN SET FORTH CUSTOMER'S SOLE AND EXCLUSIVE REMEDY IN THE EVENT OF A DEFECT IN WORKMANSHIP OR MATERIALS AND SHALL NOT BE DEEMED TO HAVE FAILED ITS ESSENTIAL PURPOSE SO LONG AS THUASNE IS WILLING AND ABLE TO CARRY OUT THE TERMS OF THIS EXCLUSIVE WARRANTY.

Customers should carefully inspect their products as soon as they receive them for manufacturer's defects. Manufacturer's defects are flaws that are existent or apparent when the package is opened. Note that manufacturer warranties do NOT include: normal wear and tear, pulled threads, runs, wear at the toe, or holes due to fingernails, toenails, or snags; or any problem caused by the wearer in application, improper use, or improper care of the garment. It is recommended that compression products be replaced every 6 months as the compression can break down with normal wear and care. These occurrences will not be considered as defects. If a customer feels that a hole, snag, tear or run is the result of a manufacturer's defect, please contact us. We take quality issues very seriously and we will do our best to resolve the matter.

Thuasne products purchased from unauthorized resellers are not covered by our warranty. Thuasne products purchased from any marketplace seller will not be covered by our warranty.



**THUASNE®**

**THERAFIRM®**

**THUASNE, LLC**

**TOLL FREE** 1.800.821.3094 (USA/Canada)

**TOLL FREE FAX** 913.281.5455 (USA)

Certified Quality System • ISO 13485

Proudly made in Hamlet, North Carolina  
and Kansas City, Kansas, USA

Contact your local Thuasne Representative

Four horizontal lines for contact information.

**ISO 13485 Certification**

is used by organizations involved in the design,  
production, installation and servicing of medical devices.



Manufactured under a BSI certified Quality Management System.  
Certification #FM584357

*View the Catalog Online*

**WHOLESALE.THERAFIRM.COM**

